 0373	Conforme alla direttiva Europea 93/42/CEE Complies with the European Directive 93/42 / EEC Cumple con la directiva europea 93/42/EEC Em conformidade com a diretiva europeia 93/42/EEC Conform directieve europene 93/42/EEC Avprva direktiv 93/42/EEC le ovomlu V souladu s evropskou směrnicí 93/42/EEC Съответств на европейската директива 93/42/EEC Megfelel a 93/42/EEC európai irányelvnek
	Sterile se non aperto o danneggiato / Raggi Beta Sterile if not opened or damaged / Beta rays Estéril si no está abierto o dañado / Beta Rays Estéril se fechado ou danificado / Beta Rays Steril dacă este nedeschis sau deteriorat/Raze Beta Ağırlımsız vəya hasar görmüşse sterildir / Beta ışınları Sterilni, paketi acılmıy şevetin nebo poškozen / Beta paprsky Стерильн, ако с неопорен или пощепен / Beta Rays Steril, ha bontatlan vagy sérült / Beta Rays
	Non riutilizzare Do not reuse No reutilizar NÃO reutilize Nu reutilizați Tekrar kullanmayın Neponožujte znovu Не використовуйте повторно Ne használja újra
	Consultare le istruzioni per l'uso Consult the instructions for use Consultar las instrucciones de uso Consulte as instruções de uso Consultati instrucțiunile de utilizare Kullannım talimatlarına bakın Kullanim için darsmanlığa danışın Prostaduje si návod k použití Olvasa el a használati utasítást
0°C 	Temperatura di conservazione Storage temperature Temperatura de almacenamiento Temperatura de depozitare Depolama sıcaklığı Skladovací teplota Temperatura na съхранение Tárolási hőmérséklet
	Numero di catalogo Catalog number Número de catálogo Número de catálogo Număr de catalog Katalóg szám Katalogové číslo Каталожен номер Katalógus szám
	Codice partita / Numero lotto Match code / lot number Código de lote / Número de lote Código do Lote / Número do Lote Cod lot / Număr lot Parti Kodu / Parti Numarası Kód sárze / Číslo sárze Код на партида / номер на партида Tételkód / téfelszám
	Data di scadenza Expiration date Fecha de caducidad Data de validade Data expirării Son kullanna tarihi Datum spotřebny Срок на годност Lejárnai dátum
	Data di fabbricazione Manufacturing date Fecha de fabricación Data de fabricaço Data fabricației Értelm tarihi Datum výroby Дата на производство Gyártás időpontja
	Non utilizzare se l'imballaggio è danneggiato Don't use if package is damaged No utilizar si el embalaje está dañado. NÃO use se a embalagem estiver danificada Nu utilizați dacă ambalajul este deteriorat Ambalaj hasarlıysa kullanmayın Neponožujte, pokud je obal poškozen Не використовуйте, ако опаківката є пошкоджена Ne használja, ha a csomagolás sérült
	Non risterrilizzare Do not reesterilize No reesterilizar NÃO reesterilizar Nu reesterilizați Yeniden sterilize etmeyin Nesterilizați Не стерилизуйте повторно Ne sterilizálja újra
	Fabricante Manufacturer Fabricante Producător Üretici firma Výrobce производитель Gyártó

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Revisione : 002 del 03/04/2023

EN Instructions for use - Read Carefully

EXAFLEX

EXAFLEX High-density collagen matrix, made up of non-reticulated bovine pericardium. It comes from cattle less than 20 months old. The product allows the formation of newly formed collagen and the subsequent vascularization of the implant. Subsequently, the implant will be remodeled and replaced with host tissue. In case of contact with viscera, the EXAFLEX device reduces adherence phenomena to a minimum.

INDICATIONS FOR USE: The EXAFLEX membrane and powder are intended for the following uses:

- urinary incontinence, pelvic floor reconstruction, vulnology
- repair of vaginal and rectal prolapse
- surgical repair of soft tissue deficiencies, placement of breast implants
- Abdominal and chest wall repair
- reinforcement of muscle flaps, ligaments and tendons
- hernia repair, vascular surgery
- use in dentistry in surgical repair of soft tissue and as a barrier against epithelial invasion of bone tissue.

CONTRAINDICATIONS: Known or presumed sensitivity to derivatives of bovine origin. The product has not been tested on pregnant subjects.

ADVERSE REACTIONS: As with any surgical procedure, adverse reactions are possible.

PRECAUTIONS: The product is disposable, therefore after opening the package, any leftovers cannot be reused or re-sterilized. Before use, check the integrity of the packaging, otherwise do not use. Improper use could cause complications of a surgical nature. Do not use the product if not stored according to the manufacturer's instructions.

ATTENTION: During its use, the product must never be dehydrated, use constant immersions in sterile physiological solution. In case you are forced to use the product in contaminated sites, it is recommended to take appropriate local and systemic measures in order to contain the contamination/infection.

STORAGE: The product must be stored in clean and dry environments. The storage temperature during the shelf life of the product must not be lower than 0°C / 32°F. DO NOT store in refrigerated cells.

TRANSPORT: Transport with intact packaging can take place under normal conditions. Both the product and the original packaging may be subjected to sudden changes in temperature for short periods.

INSTRUCTIONS: The instructions below are not guidelines for surgical technique and are not a substitute for institutional protocols or professional clinical judgments regarding patient care. Check that the product packaging is in good condition and that the packages containing the membrane are intact. Open the packaging respecting the asepsis of the operating room. Rehydrate the membrane by immersion in sterile saline at room temperature. The product can be soaked even for long periods, 2 or more hours. The product can be immersed in antibiotic-loaded physiological solution. Use the product according to surgical needs.

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In case of contact with viscera, the EXAFLEX device reduces adherence phenomena to a minimum.

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