 0373	Conforme alla direttiva Europea 93/42/CEE Complies with the European Directive 93/42 / EEC Cumple con la directiva europea 93/42/EEC Em conformidade com a diretiva europeia 93/42/EEC Conform directive europene 93/42/EEC Αvαρρ dιrεktív 93/42/EEC lie vovmlu V souladu s evropskou směrnicí 93/42/EEC Сvответств на европейската директива 93/42/EEC Megfelel a 93/42/EEC európai irányelvnek	
	Sterile se non aperto o danneggiato / Raggi Beta Sterile if not opened or damaged / Beta rays Estéril si no está abierto o dañado / Beta Rays Estéril se fechado ou danificado / Beta Rays Steril dacă este nedeschis sau deteriorat/Raze Beta Aqlınamıssa veyta basar gırımlısse sterildir / Beta ışınları Sterilni, pakoval iseni otvoren nebo polozoen / Beta papsky Стерильн, ако с неопрепен тум /стерилизиран / Beta Rays Steril, ha bontalan vagy sérült / Beta Rays	
	Non riutilizzare Do not reuse No reutilizar Năo reutilize Nu reutilizati Tekrar kullanmayın Неповъзвeжe знову Не използвайте повторно Ne használja újra	
	Consultare le istruzioni per l'uso Consult the instructions for use Consultar las instrucciones de uso Consulte as instruções de uso Consultati instrucțiunile de utilizare Kullannim talimurlarına bakın Kullanim için dānsımanlığa dānşın Prostředje si návod k použití Olvassa el a használati utasítást	
	Temperatura di conservazione Storage temperature Temperatura de almacenamiento Temperatura de armazenamıento Temperatura de depozitare Depolama sıcaklıđı Skladovaci teplota Temperatura na sklārenāne Tárolási hőmérséklet	
	Numero di catalogo Catalog number Numero de catalogo Número de catálogo Număr de catalog Katalođni numaras Katalođovo číslo Katalođovı numere Katalođovı szám	
	Codice partita / Numero lotto Match code / lot number Código de lote / Número de lote Código do Lote / Número do Lote Cod lot / Număr lot Parti Kodu / Parti Numarası Kód sárze / Culo sárze Kod na partıdla / номер на партида Tételkód / tételszám	
	Data di scadenza Expiration date Fecha de caducidad Data de validade Data expirării Son kullanna tarihi Datum spotteby Срок на годност Lejárati dátum	
	Data di fabbricazione Manufacturing date Fecha de fabricación Data de fabricađão Data fabricării Értelmı tarihi Datum výroby Дата на произвoдство Gyártás időpontja	
	Non utilizzare se l'imballaggio è danneggiato Don't use if package is damaged No utilizar si el embalaje está dañado. Năo use se a ambalagem estiver danificada Nu utilizați dacă ambalajul este deteriorat Ambalaj hasarıysa kullannmayın Неповъзвeжe, pokud je obal poškozen Не използвайте, ако опаковката е повредена Ne használja, ha a csomagolás sérült	
	Non risterrilizzare Do not resterilize No reesterilizar Năo reesterilizar Nu resterilizați Yeniden sterilize etmeyin Nesterilizație Не стерилизайте повторно Ne sterilizálja újra	
	Fabricatore Manufacturer Fabricante Producător Üretici firması	Výrobce производител Gyártó

IT

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EN

Instructions for use - Read Carefully

EXAFLEX Dust

EXAFLEX High-density collagen matrix, made up of non-reticulated bovine pericardium. It comes from cattle less than 20 months old.

The product allows the formation of newly formed collagen and the subsequent vascularization of the implant. Subsequently, the implant will be remodeled and replaced with host tissue.

In case of contact with viscera, the *EXAFLEX* device reduces adherence phenomena to a minimum.

INDICATIONS FOR USE:

EXAFLEX powder is intended for the following uses:

- urinary incontinence, pelvic floor reconstruction, vulnology
- repair of vaginal and rectal prolapse
- regeneration of the dermis after traumatic events and post-traumatic injuries
- treatment of imperfections of the dermis
- correction of blemishes such as skin hypotrophy and hypotonia
- restoration and correction of scars
- surgical repair of tunneled fistulas.

CONTRAINDICATIONS: Known or presumed sensitivity to derivatives of bovine origin. The product has not been tested on pregnant subjects. Do not use in the periorbital and glabellar area. Do not use on patients with collagen disease, Connective tissue diseases, porphyria and on children.

ADVERSE REACTIONS: As with any surgical procedure, adverse reactions are possible.

PRECAUTIONS: The product is disposable, therefore after opening the package, any leftovers cannot be reused or re-sterilized. Before use, check the integrity of the packaging, otherwise do not use. Improper use could cause complications of a surgical nature. Do not use the product if not stored according to the manufacturer's instructions.

ATTENTION: During its use, the product must never be dehydrated, use sterile physiological solution. In case you are forced to use the product in contaminated sites, it is recommended to take appropriate local and systemic measures in order to contain the contamination/infection.

STORAGE: The product must be stored in clean and dry environments. The storage temperature during the shelf life of the product must not be lower than 0°C / 32°F. DO NOT store in refrigerated cells.

TRANSPORT: Transport with intact packaging can take place under normal conditions. Both the product and the original packaging may be subjected to sudden changes in temperature for short periods.

INSTRUCTIONS: The instructions below are not guidelines for surgical technique and are not a substitute for institutional protocols or professional clinical judgments regarding patient care. Check that the product packaging is in good condition and that the packages containing the membrane are intact. Open the packaging respecting the asepsis of the operating room. The product can be rehydrated with physiological solution added with medicines or anesthetics.

Use the product according to surgical needs

EXAFLEX Preparation

Always handle the product using aseptic techniques.

Check that the package is intact including the protective cap of the syringe
Containing pericardium powder. Do not use if the packaging is damaged or the cap is missing.

Fill the empty syringe with the prescribed saline according to the product code in use.
Remove the cap of the pericardium powder syringe and screw on the Luer connector Lock to the syringe.

Attach the syringe containing saline to the other end of the Luer Lock connector.

Hold both syringes firmly in your hands.

Squeeze all of the saline into the pericardium powder syringe.

Depress the plungers back and forth at least 10 times to prepare the pericardial suspension.

Check that the product is completely mixed, then collect it totally in only one of the two syringes.

Remove the Luer Lock connector and empty syringe.

Attach a needle or silicone injector.

The silicone injector can be shortened with scissors.

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